APO A-I

# Apolipoprotein A-I

## Turbidimetry

### Quantitative determination of apolipoprotein A-I (APO A-I) IVD

Store 2 - 8°C.

#### PRINCIPLE OF THE METHOD

Turbidimetric test for the measurement of apolipoprotein A-I in human serum or plasma.

Anti- Apo A-I antibodies when mixed with samples containing Apo A-I, form insoluble complexes. These complexes cause an absorbance change, dependent upon the Apo A-I concentration of the patient sample, that can be quantified by comparison from a calibrator of known Apo A-I concentration.

#### CLINICAL SIGNIFICANCE<sup>1</sup>

Apo A-I is the major structural apolipoprotein in HDL and constitutes about 70% of the total protein. Apo A-I is a cofactor for lecithin-cholesterol-acyl-transferase (LCAT), the enzyme responsible for forming cholesteryl esters in plasma and plays an important role in the transport of cholesterol from peripheral tissues to the liver, to be finally excreted. Measurements of Apo A-I concentration is specially important in detecting coronary heart disease risk (CHD) as well as in the diagnostic of hyperlipoproteinemia. Concentrations < 120 mg/L are associated to an increased CHD risk, while concentrations ≥ 160 mg/L may even protect from the same risk. Patients with deficiencies in Apo A-I synthesis may highly increase the CHD risk. Tanger disease, a consequence of an Apo A-I catabolism defect, is characterized by several reduced plasma HDL cholesterol (HDL-c) concentration, abnormal HDL composition and accumulation of cholesteryl esters in many body tissues. Plasma HDL-c and Apo A-I concentrations in homozygotes are very low, while Apo A-I concentration is less than 10% of its normal concentration. Heterozygotes are characterized by half-normal concentration of HDL-c, Apo AI and Apo -II. Current evidence suggests that these patients have increased incidence of CHD.

#### REAGENTS

<b>Diluent (R1)</b>	Tris buffer 20 mmol/L, PEG, pH 8.3. Preservative.
<b>Antibody (R2)</b>	Goat serum, anti-human Apo B, tris 50 mmol/L, pH 7.5. Preservative.
<b>Optional</b>	APO CAL ref: 93005

#### CALIBRATION

The assay and the value of the calibrator concentration have been standardized against the Certified Reference Material WHO/IFCC SP1-01 (CDC, USA). It is recommended the use of the APO CAL Calibrator for calibration.

#### PREPARATION

**Reagents:** Ready to use.

**Calibration Curve:** Prepare the following APO CAL Calibrator dilutions in NaCl 9 g/L as diluent. Multiply the concentration of the Apo A-I calibrator by the corresponding factor stated in table bellow to obtain the Apo A-I concentration of each dilution.

Calibrator dilution	1	2	3	4	5	6
Calibrator (µL)	--	10	25	50	75	100
NaCl 9 g/L (µL)	100	90	75	50	25	-
Factor	0	0.1	0.25	0.5	0.75	1.0

#### STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date.

**Reagent deterioration:** The presence of particles and turbidity.

Do not freeze; frozen Antibody or Diluent could change the functionality of the test.

#### ADDITIONAL EQUIPMENT

- Spintech 240 autoanalyzer
- Laboratory equipment.

#### SAMPLES

Fresh serum or plasma. EDTA or heparin should be used as anticoagulant. Stable 2 weeks at 2-8°C or 3 months at -20°C.

The samples with presence of fibrin should be centrifuged before testing.

Do not use highly hemolyzed or lipemic samples.

#### QUALITY CONTROL

Control sera are recommended to monitor the performance of manual and automated assay procedures. Spinreact Apolipoprotein Control (Ref.:93006) is available. Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

#### REFERENCE VALUES<sup>5</sup>

Between 122 – 161 mg/dL.

Each laboratory should establish its own reference range.

#### SPINTECH 240 APPLICATION

Item Name APO A1			
<b>DATA INFORMATION</b>		<b>CALIBRATION</b>	
Units	mg/dL	TYPE	Spline
Decimals	0		
<b>ANALYSIS</b>		<b>STANDARD</b>	
Type	END	#1 0.10 x Cal. Val	#4 0.75 x Cal. Val
W.Length 1	340	#2 0.25x Cal. Val	#5 1.00 x Cal. Val
		#3 0.50 x Cal. Val	#6
		<b>NORMAL RANGE (37°C)</b>	
Method	Turbidimetry	LOW	HIGH
		SERUM	MALE
			FEMALE
		URINE	
SLOPE	INTER		
1.000 x +	0		
Item Name APO A1			
<b>ASPIRATION</b>		<b>DATA PROCESS</b>	
KIND	Single <input type="checkbox"/> Double <input checked="" type="checkbox"/>	<b>ABSORBANCE LIMIT</b>	
		READ	LOW -3.000
		START	END
SAMPLE	2 µL	MAIN	42 43
REAGENT 1	240 µL	SUB	30 31
REAGENT 2	60 µL		
		<b>ENDPOINT LIMIT 3</b>	
		<b>LINEAR CHECK (%)</b>	
		<b>FACTOR</b>	
Third Mix	<input checked="" type="checkbox"/> OFF <input type="checkbox"/> ON	Blank Correction 1.000	
R1 Blank	<input checked="" type="checkbox"/> Water <input type="checkbox"/> R1-B		
<b>MONITOR</b>		<b>PROZONE CHECK</b>	
0 LEVEL POINT	1	START	END
SPAN	3.000	LIMIT (%)	
		FIRST	
		SECOND	<input checked="" type="checkbox"/> Low High
		THIRD	<input checked="" type="checkbox"/> Low High

\*\* Modify reagents and sample volumes according to the range accepted but keeping always the mentioned ratio.

#### INTERFERENCES

Hemoglobin (20 g/L), bilirubin (40 mg/dL), lipemia (< 5 g/L), and rheumatoid factor (800 IU/mL) do not interfere. Other substances may interfere <sup>6,7</sup>.

#### PERFORMANCE CHARACTERISTICS

- 1. Measurement range:** Up to 250 mg/dL, under the described assay conditions. Samples with higher concentrations, should be diluted 1/5 in NaCl 9 g/L and retested again. The linearity limit depends on the sample / reagent ratio. It will be higher by decreasing the sample volume, although the sensitivity of the test will be proportionally decreased.
- 2. Detection Limit:** Values less than 0,1 mg/dL give non-reproducible results.
- 3. Precision:** The reagent has been tested for 20 days, using three levels of serum in a EP5-based study (NCCLS).

EP5	CV (%)		
	27,22mg/dL	65,74 mg/dL	131,07 mg/dL
Total	4%	3,7%	4,8%
Within Run	2,2%	0,8%	1,1%
Between Run	2,3%	1,3%	1,4%
Between Day	2,4%	3,3%	4,5%

- 4. Accuracy:** Results obtained using this reagent (y) were compared to those obtained with a Bayer immunoturbidimetric method. 39 samples ranging from 50 to 200 mg/dL of Apo A-I were assayed. The correlation coefficient (r) was 0,92 and the regression equation  $y = 1,18x - 37,8$ .

The results of the performance characteristics depend on the used analyzer.

#### NOTES

Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

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