Creatine kinase - MB

Anti CK-M. Immunoinhibition. Kinetic UV. Liquid

Quantitative determination of creatine kinase MB (CK-MB)
IVD

Store at 2-8°C

PRINCIPLE OF THE METHOD

The procedure involves measurement of CK activity in the presence of an antibody to CK-M monomer. This antibody completely inhibits the activity of CK-MM and half of the activity of CK-MB while not affecting the B subunit activity of CK-MB and CK-BB. Then it's used the CK method to quantitatively determine CK-B activity^{1,2}. The CK-MB activity is obtained by multiplying the CK-B activity by two.

CLINICAL SIGNIFICANCE

CK-MB is an enzyme formed by the association of two subunits from muscle (M) and nerve cells (B). CK-MB is usually present in serum at low concentration; it is increased after an acute infarct of myocardium and later descends at normal levels. Also is increased, rarely, in skeletal muscle damage^{5,6}.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

REAGENTS

R 1	Imidazol, pH 6.7 D-Glucose N-Acetyl-L-Cysteine Magnesium acetate NADP EDTA Hexokinase	125 mmol/L 25 mmol/L 25 mmol/L 12,5 mmol/L 2,52 mmol/L 2,02 mmol/L ≥6 800 U/L
Anti-human polyclonal CK-M antibody (sheep) sufficient to inhibit up to 2 000 U/L of CK-MM		
R 2	ADP AMP di-Adenosine-5'-pentaphosphate Glucose-6-phosphate DH Creatine phosphate	15,2 mmol/L 25 mmol/L 103 mmol/L ≥8 800 U/L 250 mmol/L

PREPARATION

All the reagents are ready to use.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented.

Do not use reagents over the expiration date.

Signs of reagent deterioration:

- Presence of particles and turbidity.
- Blank absorbance (A) at 340 nm ≥ 1.60.

ADDITIONAL EQUIPMENT

- SPIN640 / SPIN640Plus Autoanalyzer.
- General laboratory equipment.

SAMPLES

Serum free of hemolysis or heparin plasma¹: Stability 7 days at 2-8°C, protected from light.

CK-MB activity decreases a 10% after 24 hours at 4°C or 1 hour at 25°C.

REFERENCE VALUES

CK-MB > 24 U/L

CK-MB Activity × 100: 6-25 % CK-MB Activity in the sample
CK Total Activity

These values are for orientation purpose; each laboratory should establish its own reference range.

QUALITY CONTROL

CK-Nac/CK-MB specific control sera (Ref. 1002260 or Ref. 1002262) are recommended to monitor the performance of assay procedures.

If control values are found outside the defined range, check the instrument, reagents and technique for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

BARCODED REAGENTS LOAD MUST BE PRECEDED OF A SPINREACT "DATABASE" COPY INTO THE ANALYZER SOFTWARE. IT IS AVAILABLE UNDER REQUEST TO SPINREACT.

SPIN640 APPLICATION

TEST INFORMATION		REAGENT VOLUME		
Nº	**	Vol. R1	240	
Test	CKMB	Vol. R2	60	
Full Name	Creatin Kinase-MB	Vol. R3		
Standard nº	1	Vol. R4		
SAMPLE VOLUME		RESULT SETUP		
Vol. Sample Stand.	12	Decimal	1	Slope 1
Vol. Sample Increas.		Unit	U/L	Inter. 0
Vol. Sample Dec				
REACTION PARAMETERS				
Reac. Type	Fixed Time	Direction	Increase	
Pri. Wave.	340	Reagent Blank	0-0	
Sec. Wave.		React. Time	50-80	

SPIN640Plus APPLICATION

EDIT PARAMETERS			
Test	CKMB	No.	**
Full name	CKMB	Print name	CKMB
Reac. Type	Fixed Time	Direction	Increase
Pri. Wave.	340	Sec. Wave.	
Unit	U/L	Decimal	0
Reagent Blank	0-0	React. Time	56-82
Vol. Sample	12 ul	R1	240 ul
Increased		R2	60 ul
Decreased		R3	
Sample blank		R4	

The Calibration is stable until 32 days. After this period the Calibration must be performed again in order to obtain good results.

The Calibration of this parameter should be performed using the CK-Nac / CK-MB CONTROL Ref. 1002260 or Ref. 1002262 (high level).

PERFORMANCE CHARACTERISTICS

Measuring range: From detection limit of 1,9 U/L to linearity limit of 318 U/L.

If the results obtained were greater than linearity limit, dilute the sample 1/1 with NaCl 9 g/L and multiply the result by 2.

Precision:

	Intra-assay		Inter-assay	
Mean (U/L)	33,7	166,5	31,3	161,0
SD	1,00	3,76	1,19	3,47
CV (%)	2,96	2,26	3,81	2,15

Sensitivity: 1 U/L = 0,000134 (A).

Accuracy: Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x).

The results obtained were the following:

Correlation coefficient (r)²: 0,999.

Regression equation: y = 0,976 x - 0,269.

The results of the performance characteristics depend on the analyzer used.

LIMITATION OF THE PROCEDURE

- 1- The method will also measure any CK-BB isoenzyme present in serum. The activity of the isoenzyme is usually negligible; however, if a significant amount of CK-BB activity is present the CK-MB activity will be overestimated.
- 2- A macro form of BB (immunoglobulin complexed) has been observed which will be measured as B in the assay. If the measured CK-B activity exceeds 20% of the total CK activity, the presence of macro BB should be suspected.

BIBLIOGRAPHY

1. Abbot B et al. Creatinine kinase. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984: 1112-116.
2. Gerhardt W et al. Creatine kinase B-Subunit activity in serum after immunohinhibition of M-Subunit activity. Clin Chem 1979;(25/7): 1274-1280.
3. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
4. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
5. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
6. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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BIBLIOGRAPHIE

- Abbot B et al. Creatinine kinase. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984: 1112-116.
- Gerhardt W. et al. Creatine kinase B-Subunit activity in serum after immunohinhibition of M-Subunit activity. Clin Chem 1979;(25/7): 1274-1280.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AAC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AAC 2001.
- Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AAC 1999.
- Tietz NW et al. Clinical Guide to Laboratory Tests, 3rd ed AAC 1995.

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Ref: MD41254	Cont.	R1: 4 x 40 mL
		R2: 2 x 20 mL

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BIBLIOGRAFIA

- Abbot B et al. Creatinine kinase. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984: 1112-116.
- Gerhardt W. et al. Creatine kinase B-Subunit activity in serum after immunohinhibition of M-Subunit activity. Clin Chem 1979;(25/7): 1274-1280.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
- Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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Ref: MD41254

Cont.

R1: 4 x 40 mL

R2: 2 x 20 mL