
Quantitative determination of cholesterol IVD

Store at 2-8°C

PRINCIPLE OF THE METHOD

The cholesterol present in the sample originates a coloured complex, according to the following reactions:

The intensity of the color formed is proportional to the cholesterol concentration in the sample^{1,2}.

CLINICAL SIGNIFICANCE

Cholesterol is a fat-like substance called a lipid that is found in all body cells. The liver makes all of the cholesterol the body needs to form cell membranes and to make certain hormones.

The determination of serum cholesterol is one of the important tools in the diagnosis and classification of lipemia.

High blood cholesterol is one of the major risk factors for heart disease^{5,6}.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

REAGENTS

R	PIPES pH 6.9	90 mmol/L
	Phenol	26 mmol/L
	Cholesterol esterase (CHE)	1000 U/L
	Cholesterol oxidase (CHOD)	300 U/L
	Peroxidase (POD)	650 U/L
	4 – Aminophenazone (4-AP)	0,4 mmol/L

PREPARATION

The reagent is ready to use.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use. Do not use reagents over the expiration date.

Signs of reagent deterioration:

- Presence of particles and turbidity.
- Blank absorbance (A) at 505 nm ≥ 0.26.

ADDITIONAL EQUIPMENT

- SPIN 800 Autoanalyzer.
- General laboratory equipment.

SAMPLES

Serum or plasma^{1,2}: Stability of the sample 7 days at 2-8°C or freezing at -20°C will keep samples stable for 3 months.

REFERENCE VALUES

Risk evaluation^{5,6}:

Less than 200 mg/dL	Normal
200-239 mg/dL	Borderline
240 mg/dL and above	High

These values are for orientation purpose; each laboratory should establish its own reference range.

QUALITY CONTROL

Control sera and calibrators are recommended to monitor the performance of assay procedures: SPINROL H Calibrator, SPINROL H Normal and Pathologic (Ref. 1002011, 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and technique for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

BARCODED REAGENTS LOAD MUST BE PRECEDED OF A SPINREACT "DATABASE" COPY INTO THE ANALYZER SOFTWARE. IT IS AVAILABLE UNDER REQUEST TO SPINREACT.

SPIN 800 APPLICATION

EDIT PARAMETERS			
Test	CHOL	No.	**
Full name	CHOLESTEROL	Print name	CHOL
Reac. Type	End Point	Direction	Increase
Pri. Wave.	505	Sec. Wave.	
Unit	mg/dL	Decimal	0.1
Reagent Blank	3 -4	React. Time	22 - 24
Vol. Sample	3 ul	Increas.	Decrease
		Decrease	Sample blank
R1	300 ul	R2	R3
			R4
CALIBRATION			
Rule	Linear two points	Water	0
		Calibrator	*

The Calibration is stable until **36 days**. After this period the Calibration must be performed again in order to obtain good results.

PERFORMANCE CHARACTERISTICS

Measuring range: From detection limit 0,00 mg/dL to linearity limit 1000 mg/dL.

If the concentration is greater than linearity limit dilute 1/2 the sample with CINa 9 g/L and multiply the result by 2.

Precision:

	Intra-assay (n=20)		Inter-assay (n=20)	
	Mean (mg/dL)	SD	CV (%)	
Mean (mg/dL)	99	201	96	197
SD	0,83	1,41	1,75	6,41
CV (%)	0,84	0,70	1,82	3,26

Sensitivity: 1 mg/dL = 0,0019 (A).

Accuracy: Results obtained using SPINREACT reagents did not show systematic differences when compared with other commercial reagent.

The results obtained using 50 samples were the following:

Correlation coefficient (r): 0,99549.

Regression equation: y=0,911x + 2,624.

The results of the performance characteristics depend on the analyzer used.

NOTES

1. LCF (Lipid Clearing Factor) is integrated in the reagent.
2. Calibration with the aqueous standard may cause a systematic error in automatic procedures. In these cases, it is recommended to use a serum Calibrator.
3. Use clean disposable pipette tips for its dispensation.

BIBLIOGRAPHY

1. Naito H.K. Cholesterol. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1194-11206 and 437.
2. Meiatini F. et al. The 4-hydroxybenzoate/4-aminophenazone Chromogenic System. Clin Chem 1978; 24 (12): 2161-2165.
3. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
4. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
5. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
6. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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