# Methadone Metabolite / Metabolito de Metadona Multilevel Calibrator Urine / Calibrador Multinivel Orina

DoA Liquid Test

## Methadone Metabolite (EDDP) Calibrator for Urine IVD

Store at 2 - 8°C.

### INTENDED USE

The Methadone Metabolite (EDDP) Multilevel Calibrator is used for the calibration of Spinreact DoA-EIA semi-quantitative liquid test for EDDP.

### REAGENTS

5 Levels Ref.: 933040	EDDP Calibrator for Urine
Level 0	<b>Negative Human Urine. Sodium azide 0.5 g/l.</b>
Level 1	<b>EDDP Calibrator Level 1. Human Urine. Sodium azide 0.5 g/l. EDDP</b>
Level 2	<b>EDDP Calibrator Level 2. Human Urine. Sodium azide 0.5 g/l. EDDP</b>
Level 3	<b>EDDP Calibrator Level 3. Human Urine. Sodium azide 0.5 g/l. EDDP</b>
Level 4	<b>EDDP Calibrator Level 4. Human Urine. Sodium azide 0.5 g/l. EDDP</b>

Each level has the following concentrations of drug:

ng/ml	Level 0	Level 1	Level 2	Level 3	Level 4
EDDP	0	200	300	600	1000

### PREPARATION AND STABILITY

The calibrators are ready to use. No calibrator preparation is required. The calibrator is stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use.

### WARNING AND PRECAUTIONS

- The calibrators are for in vitro diagnostic use only. Harmful if swallowed.
- The calibrators contain sodium azide that may react with lead or copper plumbing to form potentially explosive metal azide. When disposing such calibrators or wastes always flush with a large volume of water to prevent azide build-up.
- Do not use the calibrators beyond their expiration dates.
- Keep all containers closed when not in use to avoid microbial contamination.
- Do not mix calibrators from different manufacturers.
- Do not freeze calibrators.

Handle all urine specimens as if they are potentially infectious.

### PROCEDURE

The calibrators should be used in the EDDP reagent for urine.

### REFERENCES

1. Urine testing for Drug of Abuse. National Institute on Drug Abuse (NIDA) Research Monograph 73, 1986.
2. European Laboratory Guidelines for Legally Defensible Workplace Drug Testing – Urine Drug Testing, EWDTs, 2002.
3. Spinreact EDDP- EIA (Refs #930078/930080) package inserts.

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