GPT (ALT)-LQ

GPT (ALT)-LQ

NADH. Kinetic UV. IFCC rec. Liquid.


Quantitative determination of alanine aminotransferase GPT (ALT) IVD

Store at 2-8°C

PRINCIPLE OF THE METHOD

Alanine aminotransferase (ALT) or Glutamate pyruvate transaminase (GPT) catalyses the reversible transfer of an amino group from alanine to α -ketoglutarate forming glutamate and pyruvate.

The pyruvate produced is reduced to lactate by lactate dehydrogenase (LDH) and NADH:

The rate of decrease in concentration of NADH, measured photometrically, is proportional to the catalytic concentration of ALT present in the sample¹.

CLINICAL SIGNIFICANCE

The ALT is a cellular enzyme, found in highest concentration in liver and kidney.

High levels are observed in hepatic disease like hepatitis, diseases of muscles and traumatism; its better application is in the diagnosis of the diseases of the liver.

When they are used in conjunction with AST aid in the diagnosis of infarcts in the myocardium, since the value of the ALT stays within the normal limits in the presence of elevated levels of AST^{1,4,5}.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

REAGENTS

R 1 Buffer	TRIS pH 7.8	100 mmol/L
	Lactate dehydrogenase (LDH)	1200 U/L
	L-Alanine	500 mmol/L
R 2 Substrate	NADH	0,18 mmol/L
	α -Ketoglutarate	15 mmol/L

PRECAUTIONS

R1: H290-May be corrosive to metals.

Follow the precautionary statements given in MSDS and label of the product.

PREPARATION

All the reagents are ready to use.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use.

Do not use reagents over the expiration date.

Signs of reagent deterioration:

- Presence of particles and turbidity.
- Blank absorbance (A) at 340 nm < 1,00.

ADDITIONAL EQUIPMENT

- Autoanalyzer Spintech 240.
- General laboratory equipment.

SAMPLES

Serum or plasma¹: Stability 7 days at 2-8°C.

Temperature conversion factors

To correct results to other temperatures multiply by:

Assay temperature	Conversion factor to		
	25°C	30°C	37°C
25°C	1,00	1,32	1,82
30°C	0,76	1,00	1,39
37°C	0,55	0,72	1,00

REFERENCE VALUES^{4,5}

	25°C	30°C	37°C
Men	up to 22 U/L	29 U/L	40 U/L
Women	up to 18 U/L	22 U/L	32 U/L

Normal newborns have been reported to show a reference range of up to double the adult, attributed to the neonate's hepatocytes. These values decline to adult levels by approximately 3 months of age.

These values are for orientation purpose; each laboratory should establish its own reference range.

APPLICATION SPINTECH 240

Item Name GPT			
DATA INFORMATION			
Units	U/L		
Decimals	0		
ANALYSIS			
Type	RATE		
W.Length 1	340		
W.Length 2			
Method	KINETIC		
CORR			
SLOPE	INTER		
1,000 x +	0		
Item Name GPT			
ASPIRATION			
KIND	Single	<input checked="" type="checkbox"/> Double	
VOLUME			
SAMPLE	25	μ L	
REAGENT 1	200	μ L	
REAGENT 2	50	μ L	
Third Mix	<input checked="" type="checkbox"/> OFF	ON	
R1 Blank	<input checked="" type="checkbox"/> Water	R1-B	
MONITOR			
O LEVEL POINT	1		
SPAN	3,000		
		DATA PROCESS	
		ABSORBANCE LIMIT	
READ		LOW	-3,000
START END		HIGH	3,000
MAIN	37 49		
SUB			
		ENDPOINT LIMIT	
		LINEAR CHECK (%) 90	
		FACTOR	
Blank Correction		1,000	
		PROZONE CHECK	
START END		LIMIT (%)	
FIRST			
SECOND		<input checked="" type="checkbox"/> Low	<input checked="" type="checkbox"/> High
THIRD		<input checked="" type="checkbox"/> Low	<input checked="" type="checkbox"/> High

Blank parameter must be performed in order to get good results in CALIB screen from main menu.

QUALITY CONTROL

Control sera are recommended to monitor the performance of assay procedures: SPINCONTROL H Normal and Pathologic (Ref. 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and technique for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

PERFORMANCE CHARACTERISTICS

Measuring range: From detection limit of 0 U/L to linearity limit of 400 U/L.

If the results obtained were greater than linearity limit, dilute the sample 1/10 with NaCl 9 g/L and multiply the result by 10.

Precision:

	Intra-assay (n=20)		Inter-assay (n=20)	
	Mean (U/L)	SD	CV (%)	
Mean (U/L)	42,0	116	41,1	115
SD	0,47	0,42	0,76	1,61
CV (%)	1,11	0,36	1,85	1,40

Sensitivity: 1 U/L = 0,00052 Δ A / min.

Accuracy: Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x).

The results obtained using 50 samples were the following:

Correlation coefficient (r)²: 0,99597.

Regression equation: $y=1,1209x + 1,390$.

The results of the performance characteristics depend on the analyzer used.

BIBLIOGRAPHY

- Murray R. Alanine aminotransferase. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1088-1090.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
- Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

PACKAGING

Ref: TK41274

Cont.

R 1: 10 x 25 mL
R 2: 10 x 7 mL

GPT (ALT)-LQ

GPT (ALT)-LQ

NADH. Kinetic UV. IFCC rec. Liquid.


Quantitative determination of alanine aminotransferase GPT (ALT) IVD

Store at 2-8°C

PRINCIPLE OF THE METHOD

Alanine aminotransferase (ALT) or Glutamate pyruvate transaminase (GPT) catalyses the reversible transfer of an amino group from alanine to α -ketoglutarate forming glutamate and pyruvate.

The pyruvate produced is reduced to lactate by lactate dehydrogenase (LDH) and NADH:

The rate of decrease in concentration of NADH, measured photometrically, is proportional to the catalytic concentration of ALT present in the sample¹.

CLINICAL SIGNIFICANCE

The ALT is a cellular enzyme, found in highest concentration in liver and kidney.

High levels are observed in hepatic disease like hepatitis, diseases of muscles and traumatism; its better application is in the diagnosis of the diseases of the liver.

When they are used in conjunction with AST aid in the diagnosis of infarcts in the myocardium, since the value of the ALT stays within the normal limits in the presence of elevated levels of AST^{1,4,5}.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

REAGENTS

R 1 Buffer	TRIS pH 7.8	100 mmol/L
	Lactate dehydrogenase (LDH)	1200 U/L
	L-Alanine	500 mmol/L
R 2 Substrate	NADH	0,18 mmol/L
	α -Ketoglutarate	15 mmol/L

PRECAUTIONS

R1: H290-May be corrosive to metals.

Follow the precautionary statements given in MSDS and label of the product.

PREPARATION

All the reagents are ready to use.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use.

Do not use reagents over the expiration date.

Signs of reagent deterioration:

- Presence of particles and turbidity.
- Blank absorbance (A) at 340 nm < 1,00.

ADDITIONAL EQUIPMENT

- Autoanalyzer Spintech 240.
- General laboratory equipment.

SAMPLES

Serum or plasma¹: Stability 7 days at 2-8°C.

Temperature conversion factors

To correct results to other temperatures multiply by:

Assay temperature	Conversion factor to		
	25°C	30°C	37°C
25°C	1,00	1,32	1,82
30°C	0,76	1,00	1,39
37°C	0,55	0,72	1,00

REFERENCE VALUES^{4,5}

	25°C	30°C	37°C
Men	up to 22 U/L	29 U/L	40 U/L
Women	up to 18 U/L	22 U/L	32 U/L

Normal newborns have been reported to show a reference range of up to double the adult, attributed to the neonate's hepatocytes. These values decline to adult levels by approximately 3 months of age.

These values are for orientation purpose; each laboratory should establish its own reference range.

APPLICATION SPINTECH 240

Item Name GPT			
DATA INFORMATION			
Units	U/L		
Decimals	0		
ANALYSIS			
Type	RATE		
W.Length 1	340		
W.Length 2			
Method	KINETIC		
CORR			
SLOPE	INTER		
1,000 x +	0		
Item Name GPT			
ASPIRATION			
KIND	Single	✓ Double	
VOLUME			
SAMPLE	25	μL	
REAGENT 1	200	μL	
REAGENT 2	50	μL	
Third Mix	✓ OFF	ON	
R1 Blank	✓ Water	R1-B	
MONITOR			
0 LEVEL POINT	1		
SPAN	3,000		
		DATA PROCESS	
		ABSORBANCE LIMIT	
READ		LOW	-3,000
START END		HIGH	3,000
MAIN	37 49		
SUB			
		ENDPOINT LIMIT	
		LINEAR CHECK (%) 90	
		FACTOR	
Blank Correction		1,000	
		PROZONE CHECK	
START END		LIMIT (%)	
FIRST		✓ Low High	
SECOND		✓ Low High	
THIRD		✓ Low High	

Blank parameter must be performed in order to get good results in CALIB screen from main menu.

QUALITY CONTROL

Control sera are recommended to monitor the performance of assay procedures: SPINCONTROL H Normal and Pathologic (Ref. 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and technique for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

PERFORMANCE CHARACTERISTICS

Measuring range: From detection limit of 0 U/L to linearity limit of 400 U/L.

If the results obtained were greater than linearity limit, dilute the sample 1/10 with NaCl 9 g/L and multiply the result by 10.

Precision:

	Intra-assay (n=20)		Inter-assay (n=20)	
	Mean (U/L)	SD	CV (%)	
Mean (U/L)	42,0	116	41,1	115
SD	0,47	0,42	0,76	1,61
CV (%)	1,11	0,36	1,85	1,40

Sensitivity: 1 U/L = 0,00052 Δ A / min.

Accuracy: Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x).

The results obtained using 50 samples were the following:

Correlation coefficient (r)²: 0,99597.

Regression equation: $y=1,1209x + 1,390$.

The results of the performance characteristics depend on the analyzer used.

BIBLIOGRAPHY

- Murray R. Alanine aminotransferase. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1088-1090.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
- Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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