HDL Cholesterol

Direct. Enzymatic colorimetric




Quantitative determination of HDL cholesterol IVD

Store at 2-8°C



PRINCIPLE OF THE METHOD

Directly determination of serum HDLc (high-density lipoprotein cholesterol) levels without the need for any pre-treatment or centrifugation of the sample^{3,5}. The assay takes place in two steps.

– 1° Elimination of lipoprotein no-HDL

– 2° Measurement of HDLc

The intensity of the color formed is proportional to the HDLc concentration in the sample.

CLINICAL SIGNIFICANCE

HDL particles are high-density lipoproteins that transport cholesterol from the body tissues to the liver. Since HDL can remove cholesterol from the arteries and carry it back to the liver for their excretion, HDL is known as "good cholesterol" because high levels are thought to lower the risk of heart disease and coronary artery disease.

A low HDL cholesterol levels, is considered a greater heart disease risk^{1,2,4}. Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

REAGENTS

R 1	N,N-bis(2-hydroxyethyl)-2-aminoethanesulphonic acid pH 6,6	100 mM
	N-(2-hydroxy-3-sulfopropyl)-3,5-dimethoxyaniline (HDAOS)	0,7 mM
	Cholesterol Esterase	≥ 800 U/L
	Cholesterol oxidase	≥ 500 U/L
	Catalase	≥ 300 U/L
R 2	Ascorbic oxidase	≥ 3000 U/L
	N,N-bis(2-hydroxyethyl)-2-aminoethanesulphonic acid pH 7,0	1,1 mmol/L
	4 – Aminoantipyrine (4-AA)	100 mM
	Peroxidase	≥ 3500 U/L

PREPARATION

- R 1 and R 2: Are ready to use.

STORAGE AND STABILITY^(Note 1)

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not freeze the reagents.

- R 1 and R 2: Once opened is stable 4 weeks at 2-8°C.

Do not use reagents over the expiration date.

Signs of reagent deterioration:

- Presence of particles and turbidity.

ADDITIONAL EQUIPMENT

- MINDRAY BS-120 / BS-200E Autoanalyzer.

- General laboratory equipment.

SAMPLES

Serum heparinized plasma or EDTA plasma. If any samples show precipitates, centrifuge before using.

Stability of the sample: 6 days at 2-8°C and 1 year when stored at -70°C.

QUALITY CONTROL

Control sera and calibrators are recommended to monitor the performance of assay procedures: SPINCONTROL H Calibrator, SPINCONTROL H Normal and Pathologic (Ref. 1002011, 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and technique for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

NOTES

1. The reagent 2 presents yellowish coloration due to the peroxidase, but it does not affect its functionality.

MINDRAY BS-120 / BS-200E APPLICATION

PARAMETERS			
Test	HDL / HDL	R1	225 / 225
Nº	**	R2	75 / 75
Full Name	HDL / HDL	Sample volume	3 / 3
Standard Nº		R1 Blank	
Reac. Type	EndPoint / Endpoint	Mixed Rgt Blank	
Pri. Wavelength	578 / 570	Linearity Range	3.0 mg/dL 120.0 mg/dL
Sec. Wavelength		Linearity Limit	*
Direction	Increase / Increase	Substrate Limit	*
Reac. Time	-1 _ 18 / -1 _ 18	Factor	*
Incuba. Time		Prozone check	*
Units	mg/dL / mg/dL	q1	q2
Precision	0.1 / 0.1	q3	q4
		PC	Abs
CALIBRATION (Cal + Rgt Blk)			
Rule	One-point Linear / Two-point Linear		
Sensitivity	1 / 1		
Replicate	2 / 2		
Interval (days)	0 / 0		
Difference Limit			
SD			
Blank Response			
Error Limit			
Correlation Coefficient			

Blank parameter must be performed in order to get good results in CALIB screen from main menu. The blank calibration is stable until **35 days**. After this period the blank parameter must be performed again in order to validate the calibration.

REFERENCE VALUES

	Men	Women
Low risk	> 50 mg/dL	> 60 mg/dL
Normal risk	35 – 50 mg/dL	45 – 60 mg/dL
High risk	< 35 mg/dL	< 45 mg/dL

These values are for orientation purpose; each laboratory should establish its own reference range.

PERFORMANCE CHARACTERISTICS

Measuring range: From *detection limit* of 9,7 mg/dL to *linearity limit* of 151 mg/dL. If the results obtained were greater than linearity limit, dilute the sample 1/2 with NaCl 9 g/L and multiply the result by 2.

Precision:

Mean (mg/dL)	Intra-assay (n=30)		Inter-assay (n=30)	
	52,4	61,6	50,1	61,6
SD	1,31	1,35	2,77	2,62
CV (%)	2,52	2,18	5,57	4,27

Sensibility: 1mg/L = 0,001897 (A).

Accuracy: Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x).

The results obtained using 54 samples were the following:

Correlation coefficient (r)²: 0,994.

Regression equation: y = 0,93x + 0,033.

The results of the performance characteristics depend on the analyzer used.

BIBLIOGRAPHY

- National Institutes of Health Consensus Development Conference Statement: Triglyceride, High Density Lipoprotein and Coronary Heart Disease. Washington D.C. Feb 26-28, 1992.
- Izawa S., Okada M., Matsui H., and Horita Y. J. Medicine and Pharmaceutical Sci., 1385 - 1388, 37 (1997).
- Shih WJ, Bachorik PS, Haga JA, Myers GL, Stein EA; Clinical Chemistry, 2000; 46:3:351 – 364
- Third Report of the National Cholesterol Education Programme (NCEP) Expert Panel on Detection, Evaluation and treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III). JAMA Publication, Vol 285, No. 19, P2486 - 2497; 2001.
- Jacobs, D. et al. In Laboratory and Test Handbook; Jacobs, D.S.; Kasten, B.L., De Mott, W.R., Wolfson, W.L., Eds; Lexi - Comp Inc: Hudson (Cleveland), 1990; P. 219.

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- National Institutes of Health Consensus Development Conference Statement: Triglyceride, High Density Lipoprotein and Coronary Heart Disease. Washington D.C. Feb 26-28, 1992.
- Izawa S., Okada M., Matsui H., and Horita Y. J. Medicine and Pharmaceutical Sci., 1385 - 1388, 37 (1997).
- Shih W.J, Bachorik PS, Haga JA, Myers GL, Stein EA; Clinical Chemistry, 2000; 46:3:351 – 364
- Third Report of the National Cholesterol Education Programme (NCEP) Expert Panel on Detection, Evaluation and treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III). JAMA Publication, Vol 285, No. 19, P2486 - 2497; 2001.
- Jacobs, D. et al. In Laboratory and Test Handbook; Jacobs, D.S; Kasten, B.L., De Mott, W.R., Wolfson, W.L., Eds; Lexi - Comp Inc: Hudson (Cleveland), 1990; P. 219.

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BIBLIOGRAPHIE

- National Institutes of Health Consensus Development Conference Statement: Triglyceride, High Density Lipoprotein and Coronary Heart Disease. Washington D.C. Feb 26-28, 1992.
- Izawa S., Okada M., Matsui H., and Horita Y. J. Medicine and Pharmaceutical Sci., 1385 - 1388, 37 (1997).
- Shih WJ, Bachorik PS, Haga JA, Myers GL, Stein EA; Clinical Chemistry, 2000; 46:3:351 - 364
- Third Report of the National Cholesterol Education Programme (NCEP) Expert Panel on Detection, Evaluation and treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III). JAMA Publication, Vol 285, No. 19, P2486 - 2497; 2001.
- Jacobs, D. et al. In Laboratory and Test Handbook; Jacobs, D.S; Kasten, B.L., De Mott, W.R., Wolfson, W.L., Eds; Lexi - Comp Inc: Hudson (Cleveland), 1990; P. 219.

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Ref:MI1001096	Cont.	R1: 4 x 30 mL
		R2: 2 x 20 mL

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