
Quantitative determination of uric acid**IVD**

Store at 2-8°C

PRINCIPLE OF THE METHOD

Uric acid is oxidized by uricase to allantoin and hydrogen peroxide ($2\text{H}_2\text{O}_2$), which under the influence of POD, 4-aminophenazone (4-AP) and 2-4 Dichlorophenol sulfonate (DCPS) forms a red quinoneimine compound:

The intensity of the red color formed is proportional to the uric acid concentration in the sample^{1,2}.

CLINICAL SIGNIFICANCE

Uric acid and its salts are end products of the purine metabolism. With progressive renal insufficiency, there is retention in blood of urea, creatinine and uric acid. Elevate uric acid level may be indicative of renal insufficiency and is commonly associated with gout^{1,5,6}. Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

REAGENTS

R 1 Buffer	Phosphate pH 7,4 2-4 Dichlorophenol sulfonate (DCPS)	50 mmol/L 4 mmol/L
R 2 Enzymes	Uricase Peroxidase (POD) Ascorbate oxidase 4 – Aminophenazone (4-AP)	60 U/L 660 U/L 200 U/L 1 mmol/L

PREPARATION

All the reagents are ready to use.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use.

Do not use reagents over the expiration date.

Signs of reagent deterioration:

- Presence of particles and turbidity.
- Blank absorbance (A) at 520 nm ≥ 0,16.

ADDITIONAL EQUIPMENT

- MINDRAY BS-120 / BS-200E Autoanalyzer.
- General laboratory equipment.

SAMPLES

- Serum or plasma¹: Stability 3-5 days at 2-8°C or 6 months at -20°C.
- Urine (24 h)¹: Stability 4 days at 15-25°C, pH >8. Dilute sample 1/50 in distilled water. Mix. Multiply results by 50 (dilution factor); If urine is cloudy; warm the specimen to 60°C for 10 min to dissolve precipitated urates and uric acid. Do not refrigerate.

REFERENCE VALUES¹

Serum or plasma:

Women	2,5 – 6,8 mg/dL	≥ 149 – 405 µmol/L
Men	3,6 – 7,7 mg/dL	≥ 214 – 458 µmol/L
Urine:	250 - 750 mg/24 h	≥ 1,49 – 4,5 mmol/24 h

These values are for orientation purpose; each laboratory should establish its own reference range.

QUALITY CONTROL

Control sera and calibrators are recommended to monitor the performance of assay procedures: SPINTROL H Calibrator, SPINTROL H Normal and Pathologic (Ref. 1002011, 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and technique for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

MINDRAY BS-120 / BS-200E APPLICATION

PARAMETERS				
Test	URI / URI	R1	180 / 180	
Nº	**	R2	180 / 180	
Full Name	URIC ACID	Sample volume	9 / 9	
Standard Nº		R1 Blank		
Reac. Type	Endpoint	Mixed Rgt Blank		
Pri. Wavelength	510 / 505	Linearity Range	0,03 mg/dL	25,00 mg/dL
Sec. Wavelength		Linearity Limit	*	
Direction	Increase	Substrate Limit	*	
Reac. Time	0_17 / -5_17	Factor	*	
Incuba. Time		Prozone check	*	
Units	mg/dL / mg/dL	q1	q2	
Precision	0,01 / 0,01	q3	q4	
		PC	Abs	
<u>CALIBRATION (Cal + Rgt Blk)</u>				
Rule	One-point Linear / Two-point Linear			
Sensitivity	1 / 1			
Replicates	2 / 2			
Interval (days)	0 / 0			
Difference Limit				
SD				
Blank Response				
Error Limit				
Correlation Coefficient				

Blank parameter must be performed in order to get good results in CALIB screen from main menu. The blank calibration is stable until 35 days. After this period the blank parameter must be performed again in order to validate the calibration.

PERFORMANCE CHARACTERISTICS

Measuring range: From detection limit of 0,01647 mg/dL to linearity limit of 40 mg/dL.

If the results obtained were greater than linearity limit, dilute the sample 1/2 with NaCl 9 g/L and multiply the result by 2.

Precision:

	Intra-assay (n=20)	Inter-assay (n=20)
Mean (mg/L)	4,46	10,37
SD	0,02	0,05
CV (%)	0,46	0,44

Sensitivity: 1 mg/dL = 0,0323 (A).

Accuracy: Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x).

The results obtained using 50 samples were the following:

Correlation coefficient (r)²: 0,99734.

Regression equation: y=0,816x + 0,319.

The results of the performance characteristics depend on the analyzer used.

NOTES

1. Calibration with the aqueous standard may cause a systematic error in automatic procedures. In these cases, it is recommended to use a serum Calibrator.
2. Use clean disposable pipette tips for its dispensation.

BIBLIOGRAPHY

1. Schultz A. Uric acid. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1261-1266 and 418.
2. Fossati P et al. Clin Chem 1980;26:227-231.
3. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
4. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
5. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
6. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995

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