β₂-m TURBI

β₂-m turbilatex

Latex turbidimetry

Quantitative determination of β₂-microglobulin (β₂-m) IVD

Store 2 - 8°C.

PRINCIPLE OF THE METHOD

The β₂-m Turbilatex is a quantitative turbidimetric test for the measurement of β₂-microglobulin (β₂-m) in human serum, plasma or urine. Latex particles coated with anti-human β₂-m are agglutinated when mixed with samples containing β₂-m. The agglutination causes an absorbance change, dependent upon the β₂-m contents of the patient sample that can be quantified by comparison from a calibrator of known concentration.

CLINICAL SIGNIFICANCE

β₂-m is a protein located on the surface of human lymphocytes and other nucleated cells. Free β₂-m is filtered by the glomerulus and subsequently reabsorbed in the proximal tubular cells. Increased urinary excretion of β₂-m is a sensitive indicator of renal insufficiency. Also, the β₂-m level in serum is a useful marker of other diseases including carcinomas, lymphoid tumors, rheumatoid arthritis and AIDS.

REAGENTS

β ₂ -m Diluent (R1)	Tris buffer 20 mmol/L, pH 8,2. Preservative.
β ₂ -m Latex (R2)	Particles coated with goat IgG anti-human β ₂ -m, pH 8,0. Preservative
β ₂ -m CAL	Calibrator. β ₂ -m concentration is stated on the vial label.

PRECAUTIONS

Components from human origin have been tested and found to be negative for the presence of HBsAg, HCV, and antibody to HIV (1/2). However, handle cautiously as potentially infectious.

CALIBRATION

Use β₂- Microglobulin Calibrator Reference 1107032.

The sensitivity of the assay and the target value of the calibrator have been standardized against the 1st International β₂-m Standard from WHO. Recalibrate when control results are out of specified tolerances, when using different lot of reagent and when the instrument is adjusted.

PREPARATION

Ready for use.

β₂-m Calibrator:

Serum method: Reconstitute (→) with 1,0 mL of distilled water. Mix gently and bring to room temperature for about 10 minutes before use.
Urine method: Dilute reconstituted calibrator 1/6 with NaCl 9 g/L (50 μL calibrator + 250 μL NaCl 9 g/L).

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during use. Do not use reagents over the expiration date.

Do not freeze; frozen Latex or Diluent could change the functionality of the test.

Reagent deterioration: Presence of particles and turbidity.

β₂-m Calibrator: Stable for 1 month at 2-8°C or 3 months at -20°C.

ADDITIONAL EQUIPMENT

- SPIN640 / SPIN640Plus autoanalyzer
- Laboratory equipment.

SAMPLES

Fresh serum. Stable 7 days at 2-8°C or 3 months at -20°C.

Fresh urine. Adjust samples to pH 7-8 by the addition of K₂HPO₄. Stable 2 days at 2-8°C or 2 months at -20°C.

The samples with particles or fibrin should be centrifuged before testing. Do not use hemolized or lipemic samples.

REFERENCE VALUES

Serum: from 1,0 to 3,0 mg/L.

Urine: from 0,1 to 0,3 mg/L.

Each laboratory should establish its own reference range

QUALITY CONTROL

Control Sera are recommended to monitor the performance of manual and automated assay procedures. It should be used the SPINREACT β₂-m Control (Ref: 1107034). Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

BARCODED REAGENTS LOAD MUST BE PRECEDED OF A SPINREACT "DATABASE" COPY INTO THE ANALYZER SOFTWARE. IT IS AVAILABLE UNDER REQUEST TO SPINREACT.

SPIN640 APPLICATION

TEST INFORMATION		REAGENT VOLUME	
No	**	Vol. R1	240
Test	B2-m	Vol. R2	60
Full Name	B2-microglobulin	Vol. R3	
Standard nº	1	Vol. R4	
SAMPLE VOLUME		RESULT SETUP	
Vol. Sample Stand.	3	Decimal	0.1 Slope 1
Vol. Sample Increas.		Unit	mg/L Inter. 0
Vol. Sample Dec			
REACTION PARAMETERS			
Reac. Type	Fixed Time	Direction	Increase
Pri. Wave.	546	Reagent Blank	14-15
Sec. Wave.		React. Time	44-65

SPIN640Plus APPLICATION

EDIT PARAMETERS			
Test	B2-M	No.	**
Full name	B2-M	Print name	B2-M
Reac. Type	Fixed Time	Direction	Increase
Pri. Wave.	546	Sec. Wave.	
Unit	mg/L	Decimal	0.1
Reagent Blank	14 -15	React. Time	51 - 72
Vol. Sample	3 ul	R1	240 ul
Increased		R2	60 ul
Decreased		R3	
Sample blank		R4	

The Calibration is stable until 31 days. After this period the Calibration must be performed again in order to obtain good results.

PERFORMANCE CHARACTERISTICS

- Linearity limit:** Up to 18 mg/L (serum) and 3 mg/L (urine), under the described assay conditions. Samples higher results should be diluted 1/5 in NaCl 9 g/L and retested again. The linearity depends on the sample-reagent ratio, as well as the analyzer used. It will be higher by decreasing sample volume, although the sensitivity of the test will be proportionally decreased.
- Detection limit:** Values less than 0,09 mg/L (serum) and 0,01 mg/L (urine) give non-reproducible results.
- Prozone effect:** No prozone effect was detected up to 100 mg/L (serum) and 20 mg/L (urine).
- Sensitivity:** Δ 0,048 A. mg/L (serum) and Δ 0,228 A. mg/L (urine).
- Precision:** The reagent has been tested for 20 days, using three different β₂-m concentrations in a EP5-based study.

EP5	CV (%)		
	+/- 1 mg/L	+/- 3,2 mg/L	+/- 8,5 mg/L
Total	4,0%	3,4%	1,7%
Within Run	2,8%	2,0%	1,2%
Between Run	1,7%	1,5%	1,2%
Between Day	2,2%	2,4%	0,0%

- Accuracy:** Results obtained using this reagent (y) were compared to those obtained using a commercial reagent (x) with similar characteristics. 36 samples of different concentrations of β₂-m were assayed. The correlation coefficient (r) was 0,97 and the regression equation y = 1,709x - 2,627. The results of the performance characteristics depend on the analyzer used.

BIBLIOGRAPHY

- Bhalla, R.B. et al. Clinical Chemistry 1983; 29: 1560.
- Malaguarnera M et al. Digestive Diseases and Sciences 1997; 42: 762-766.
- Chironna et al. Int J Clin Lab Rws 1994; 24: 90-93.
- Wibell L et al. Nephron 1973; 10: 320-331.
- Berggard B et al. Scand J Clin Lab Invest 1980; 40: 13-25.
- Davey P G et al. Clin Chem 1982; 28/6: 1330-1333.
- Young DS. Effects of drugs on clinical laboratory test, 4th ed. AACC Pres, 1995.

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EDIT PARAMETERS			
Test	B2-M	No.	**
Full name	B2-M	Print name	B2-M
Reac. Type	Fixed Time	Direction	Increase
Pri. Wave.	546	Sec. Wave.	
Unit	mg/L	Decimal	0.1
Reagent Blank	14-15	React. Time	51-72
Vol. Sample	3 ul	R1	240 ul
Increased		R2	60 ul
Decreased		R3	
Sample blank		R4	

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BIBLIOGRAFIA

- Bhalla, R.B. et al. Clinical Chemistry 1983; 29: 1560.
- Malaguarnera M et al. Digestive Diseases and Sciences 1997; 42: 762-766
- Chironna et al. Int J Clin Lab Rws 1994; 24: 90-93
- Wibell L et al. Nephron 1973; 10: 320-331
- Berggard B et al. Scand J Clin Lab Invest 1980; 40: 13-25.
- Davey P G et al. Clin Chem 1982; 28/6: 1330-1333
- Young DS. Effects of drugs on clinical laboratory test, 4th ed. AACCC Press, 1995.

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Ref.:MD1107130	Cont.	R1. Diluente: 2 x 40 mL
		R2. Látex: 1 x 20 mL
		β ₂ -m -CAL 1 x 1 mL